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
Eumed S.r.l.
Via Tommaso, 68 - 35131 Padova (ITALIA)

EN

Actycoll - membranes and fleeces based on Type I collagen of equine origin for guided tissue regeneration

INSTRUCTIONS FOR USE

Actycoll - Sterile, single-use device

Actycoll Fleece
Actycoll Membrana
Actycoll RS
Actycoll GCR

INTENDED USE

Actycoll is a line of fleeces and membranes based on Type I, native, lyophilized, equine collagen for the guided regeneration of connective tissues in dental and oral-maxillofacial surgery, orthopaedics, traumatology and neurosurgery.

DESCRIPTION

Collagen is the main protein in the connective tissues of animals, representing about 25% of the total protein mass in mammals. The structural unit of collagen is composed of tropocollagen, a protein formed by three left-handed polypeptide chains that intertwine to form a triple right-handed helix. Type I collagen, consisting of two alpha 1 and one alpha 2 chains, constitutes 90% of total collagen and makes up the main connective tissues, such as skin, tendons, ligaments, bones and cornea. The specific porous structure and fibrillar nature of Actycoll fleeces and membranes guarantee adequate absorption properties of biological fluids and blood once implanted in the graft site, offering optimal support and protection during the reparative processes of connective tissues.

AVAILABLE FORMATS

The devices of the Actycoll line are available in various formats, in single and multiple packs, in the form of fleeces (Actycoll Fleece) and membranes (Actycoll Membrana, Actycoll RS and Actycoll GCR).

PRODUCT CONSTITUENTS

Actycoll Fleece
Type I equine collagen fleece.
Actycoll Membrana
Type I equine collagen membrane, smooth surfaces.
Actycoll RS / Actycoll GCR
Type I equine collagen membrane, dual-side, with a smooth and rough surface.

INDICATIONS

Actycoll is a line of membranes and fleeces based on Type I collagen of equine origin for guided tissue regeneration, indicated for the containment, protection and stabilization of bone grafts in guided bone regeneration operations or as a membrane for tissue regeneration in cartilage regeneration surgeries.

Actycoll Fleece

Fleece, capable of effectively absorbing biological fluids, is used as a haemostatic biomaterial or carrier to maintain cellular concentrates and/or autologous platelets in the implant site. Once cut to size, it can be used to fill post-extractive alveoli. Fragmented and associated with bone grafts may also act as a graft-extender to increase the volume of the material to be grafted.

Actycoll Membrana / Actycoll RS

Membranes positioned to cover bone grafts in Guided Bone Regeneration (GBR) procedures to protect the bone regeneration site from infiltration of epithelial or fibrous tissue. Actycoll membrane is used to protect post-extractive socket to reduce bone resorption, to protect the Schneider membrane in case of small ruptures or in maxillary sinus elevations. Actycoll RS, thanks to the presence of a rough side to be placed in contact with the bone defect, allows for a more natural stabilization of bone grafts.

Actycoll GCR

Membrane indicated in guided tissue regeneration procedures, for the protection of chondral defects in the treatment of grade III-IV focal joint cartilage lesions according to Outerbridge classification, with dimensions between 2 and 8 cm², after preparation of the focal cartilage lesion with micro-fractures according to Steadman's technique or nano-fractures. The membrane, positioned with the rough side facing the cartilaginous lesion, helps the stabilization of the blood clot, prevents the intra-articular dispersion of autologous cells from the bone marrow, acting as a guide for the healing of the lesion.

CONTRAINDICATIONS AND SIDE EFFECTS

Actycoll is a biocompatible device. Latex and Phthalates free. Actycoll should not be used in patients with proven individual hypersensitivity to equine collagen. It is also contraindicated in the presence of inflammation or infection at the implant site or in tissues adjacent to it.

INSTRUCTIONS FOR USE

The following instructions in no way replace the usual surgical procedures in the treatment of bone and cartilage defects.

Actycoll Fleece

The fleece used as a haemostatic can be cut out as needed and positioned directly on the implant site, without any rehydration. Even when used as a support for cellular concentrates and/or autologous platelets, the fleece should not be rehydrated, but directly combined with the selected biological fluid. Once soaked, wait 1-2 minutes and proceed with its grafting. If used as a graft-extender, Actycoll Fleece can be minced into small fragments, rehydrated with saline or autologous blood, squeezed to remove excess fluid, and then mixed with the bone grafts before implantation.

Actycoll Membrane / Actycoll RS

Once shaped as needed, Actycoll Membrane and Actycoll RS, handled with sterile surgical forceps, can be immersed in saline solution for a few seconds and directly placed in the implant site. The immersion time can be extended by a few more seconds if further softening of the membranes is desired, making them more easily adaptable to the specific conformation of the implant site. In case of abundant bleeding at the implant site, the membranes can be placed without rehydration. Actycoll Membrane and Actycoll RS do not require fixation in order to be stabilized, thanks to the natural adhesiveness of collagen-based materials. Actycoll RS should be placed with the rough side facing the bone defect in order to further facilitate stabilization.

Actycoll GCR

Application of the membrane by arthrotomy or mini-arthrotomy: clean the lesion by removing the degenerated cartilaginous tissue with a curette until the healthy cartilage is reached, taking care to create vertical and well-defined wound edges. Perform micro-fractures according to Steadman's technique, or nano-fractures, in order to promote bleeding of the subchondral bone. Cut the membrane according to the size of the defect, leaving it slightly undersized so that, once positioned with the rough side facing the bone surface, it enters completely inside the defect. It is recommended to use fibrin glue along the perimeter of the membrane to seal the edges, improving its stability.

Application of the membrane by arthroscopy: clean the lesion by removing the degenerated cartilaginous tissue with a curette until the healthy cartilage is reached, taking care to create vertical and well-defined wound edges. Cut the membrane according to the size of the defect, leaving it slightly undersized. For the subsequent surgical steps it is recommended to close the joint irrigation and insufflate CO₂ so as to operate in a dry environment while maintaining an excellent arthroscopic view. Then proceed with micro-drilling of the exposed subchondral bone surface or perform micro-fractures according to Steadman's technique. Thread the membrane through the arthroscopic access and place it onto the defect, with the rough side facing the bone surface of the defect. Spread the membrane gently with a feeler, taking care to position it completely inside the edges of the defect. Seal the edges of the defect by injecting fibrin glue along the perimeter of the membrane and wait about 2 minutes for it to stabilize after polymerization. Stop the flow of CO₂ and reactivate the joint irrigation. Perform gentle flex-extension movements to verify the membrane is tight at the lesion site.

Degradation and reabsorption, in the absence of external factors that contribute to accelerate or delay the time, usually occurs between 4-6 weeks for Actycoll Fleece, Membrane and Membrane RS, and 10-12 weeks for Actycoll GCR.

WARNINGS AND PRECAUTIONS

Read the instructions carefully before use. Actycoll is a Class III, single-use, single-patient medical device; it cannot be reused or re-sterilized. Any unused product residue must not be stored and used in other operations or for other patients and must be disposed of as contaminated hospital waste in accordance with hospital regulations.

In the event that the disposable device is mistakenly reused, as well as affecting its performance, possible infections could occur to the patient due to the lost sterility and/or possible contamination of the device itself. Failure to comply with this indication will result in a different use than that intended by the manufacturer, who can no longer guarantee the conditions of hygiene, functionality, and safety once the packaging has been opened.

Actycoll has not been tested on pregnant patients. Do not use the device if the package is open or damaged. Do not use Actycoll after the expiration date. The expiration date refers to the product in its undamaged and regularly stored packaging. Keep out of the reach of children.

PRODUCT PACKAGING

Actycoll is supplied STERILE.
Actycoll Fleece
Actycoll Membrana
Actycoll RS
Actycoll GCR

Laminated polyester-aluminium (ALU/OPA-PE) envelope, single or double. Illustrative leaflet. For all formats: the outer label of the product is placed on the outer box (final packaging), while on the outside of each envelope contained in the box there are 4 detachable adhesive patient labels made for traceability, to be applied on the patient record or other documentation.
Wrapping breakage and disposal of the package

Do not use the product if the packaging is damaged. The packaging materials do not require special disposal conditions.

STERILIZATION

Actycoll medical devices are sterilized by using 25 kGy gamma rays.

STORAGE

The devices must be stored in a cool and dry place at a temperature between 4 and 40°C avoiding direct exposure to sunlight. If stored properly and in intact packaging, the sterility of the devices is guaranteed for 5 years from the date of production (see expiration date on the outer label).

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Tiss’You S.r.l.