Guided-Drill-Stop-Control

Guided Drill System for Safe Implant Site Preparation | developed with PD Dr. Dr. Florian Bauer

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The Guided-Drill-Stop system is used for the easy and safe preparation of the implant site. The special feature of the system are the specially developed twist drills with guiding tip, with which an intermediate step is added prior to the actual extension drilling. With this step, the upper part of the implant site is expanded to the desired diameter, while the lower part is used to guide the drill. This avoids deviation from the desired axis. In addition, the stop sleeves provide the defined depth stop during all drillings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

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
Führungshülse
Drill sleeve

Surgical Guide Fixation Screw
Surgical Guide Fixation Screw

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Screw and sleeve for safe intraoral fixation of surgical guides, screw self-drilling and self-tapping, manual insertion or by means of screw driver (e.g. TL0T1).

Fig.	DBH06	GF080
L ₁ (mm)	6.0	23.0
L ₂ (mm)	5.5	6.2
□ (mm)	3.8	2.0
□ (mm)	2.2	-

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GBMN
MEMBER OF THE GLOBAL
BONE MANAGEMENT® &
RESTORATIVE NETWORK

Developed with
PD Dr. Dr. Florian Bauer
Miesbach, Germany

Praxisklinik MKG Miesbach
www.mkg-miesbach.de

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© PD Dr. Dr. Florian Bauer