



# RTRO Kit NEW

Root Tip Resection and Osteotomy | developed with PD Dr. Dr. Florian Bauer

Das RTRO Kit, welches von PD Dr. Dr. Florian Bauer zusammengestellt wurde, enthält eine abgestimmte Auswahl an Diamantinstrumenten und Hartmetallfräsern zur Durchführung einer Wurzelspitzenresektion und Osteotomie. Mit den Instrumenten HM141A und 801 wird zuerst die Wurzelspitze freigelegt, indem das darüberliegende Weichgewebe und der Knochen abgetragen werden. Die Entfernung der Wurzelspitze und des umliegenden Gewebes erfolgt anschließend mit den Instrumenten 837L, 850G, HM165 und HM166A.

The RTRO Kit by PD Dr. Dr. Florian Bauer contains a sophisticated selection of diamond instruments and carbide burs for performing root tip resection (apicoectomy and osteotomy). The HM141A and 801 instruments are used to first expose the root apex by removing the overlying bone tissue. Whereas the 837L, 850G, HM165 and HM166A instruments are applied in the subsequent removal of the root tip as well as of inflamed tissue in its vicinity.



Art.-No. RTR00

	Hartmetallinstrumente Tungsten carbide burs				Diamantinstrumente Diamond instruments			
Fig.	HM165**	HM166A**	HM141A**	HM141A**	837L	850G	801*	801
Shank <sup>1</sup>	104	104	104	104	104	104	104	104
Size <sup>2</sup>	023	021	023	035	016	018	023	033
Length mm	7.0	10.5	2.1	3.3	8.0	8.0	2.0	3.0
⌀ mm	2.3	2.1	2.3	3.5	1.6	1.8	2.3	3.3

<sup>1</sup> 104=HP

<sup>2</sup> Größter Arbeitsteildurchmesser in 1/10 mm <sup>3</sup> Largest working part diameter in 1/10 mm

⌀ Außendurchmesser External diameter



## AUF EINEN BLICK

- Abgestimmte Zusammenstellung von Diamantinstrumenten und Hartmetallfräsern für die Wurzelspitzenresektion und Osteotomie
- Alle rotierenden Instrumente in einem Kit für ein maximal effizientes Arbeiten im klinischen Alltag
- Sicheres Freilegen und Entfernen der Wurzelspitze durch hochpräzise Instrumente

## AT A GLANCE

- Sophisticated combination of high-quality diamond and carbide burs for apicoectomy and osteotomy
- All rotary instruments for exposure and tissue removal in one single kit
- Safe and efficient exposure and removal of the root tip by high-precision instruments



Praxisklinik MKG Miesbach  
www.mkg-miesbach.de



MEMBER OF THE GLOBAL  
BONE MANAGEMENT® &  
RESTORATIVE NETWORK

Developed with  
**PD Dr. Dr. Florian Bauer**  
Miesbach, Germany

