


# Close to you.


MANUFACTURER **SINCE 1921**

## Kirsch Sinus Lift Instrument Set SMALL

for sinus floor elevation, modif. acc. to Dr. Kai Zwanzig



**ZEPF** | *dental*  
newsticker



**24.995.00S**

Kirsch Sinus Lift Instrument Set SMALL,  
modif. acc. to Dr. Kai Zwanzig,  
incl. 1 x washbasket with lid 85.195.00 and  
arrangement for Sinus Lift Instrument Set

Complete Set consisting of:



**85.251.04**

Mixing Cup, stainless steel,  
with plastic lid, Ø 40 mm

**85.251.14**

Mixing Cup without plastic lid,  
stainless steel, Ø 40 mm



**41.868.01S**

**Kirsch**, Sinus 1 SMALL, Elevator 25°,  
Titanium **ZEPF**-Line Handle, double-ended



**41.868.02S**

**Kirsch**, Sinus 2 SMALL, Elevator 60°/60°,  
Titanium **ZEPF**-Line Handle, double-ended



**41.868.03S**

**Kirsch**, Sinus 3 SMALL, Elevator 90°,  
Titanium **ZEPF**-Line Handle, double-ended



**41.868.04S**

**Kirsch**, Sinus 4 SMALL, Elevator 120°/120°,  
Titanium **ZEPF**-Line Handle, double-ended

## Kirsch Sinus Lift Instrument Set SMALL

for sinus floor elevation, modif. acc. to Dr. Kai Zwanzig

Responding to numerous customer requests, the proven Kirsch shapes were minimized during the further development of the new sinus lift instruments SMALL.

- The smaller instruments allow the preparation of a minimized window.
- Regarding the angles of working tips, the previous Kirsch Instrument Set left nothing to be desired.

- The minimized analogous instruments now enable the practitioner to perform the preparation in a smaller area without injuring the Schneider's membrane.
- The **ZEPF**-Line-handle allows a safe, tactile guidance of the instruments.

Optionally, two instruments 41.868.05SL and 41.868.06SL with a longer 1st shaft are offered.



**41.868.05S** Kirsch, Sinus 5 SMALL, Elevator 90°/90°/60°;  
Titanium **ZEPF**-Line Handle, double-ended



**41.868.06S** Kirsch, Sinus 6 SMALL, Elevator 90°/90°;  
Titanium **ZEPF**-Line Handle, double-ended



**41.868.07S** Kirsch, Sinus 7 SMALL, Plugger Ø 5.0 mm, Spoon 8 x 10 mm,  
Titanium **ZEPF**-Line Handle, double-ended



**41.868.08S** Kirsch, Sinus 8 SMALL, convex Plugger Ø 4.5 mm, plane Plugger Ø 2.7 mm,  
Titanium **ZEPF**-Line Handle, double-ended



**41.747.03S** Universal Raspatory, Sinus 0 SMALL, exchangeable inserts



## Kirsch Sinus Lift Instrument Set SMALL for sinus floor elevation, modif. acc. to Dr. Kai Zwanzig

### Optional Instrumentary



**41.868.05SL**

**Kirsch**, Sinus 5 SMALL, long shank,  
Elevator 90°/90°/60°,  
Titanium **ZEPF**-Line Handle, double-ended

**41.868.06SL**

**Kirsch**, Sinus 6 SMALL, long shank,  
Elevator 90°/90°,  
Titanium **ZEPF**-Line Handle, double-ended

