

# Cytoplast™ RTM Collagen

Type I bovine collagen membrane



15 mm x 20 mm

| RTM1520 (2 per box)

actual size



20 mm x 30 mm

| RTM2030 (2 per box)



30 mm x 40 mm

| RTM3040 (2 per box)



“...I am impressed with its **handling**, but most importantly, I am impressed with its **results**.”

*Jerald Rosenberg, DMD  
Periodontist*

## MANUFACTURED FROM HIGHLY PURIFIED TYPE 1 BOVINE ACHILLES TENDON

Seguro para o paciente.

## 26 – 38 WEEK RESORPTION TIME

O longo período de reabsorção previsível limita o risco de perda de partículas por reabsorção prematura.

## HIGH TENSILE STRENGTH


Pode suturar ou prender a membrana no lugar sem rasgar.

## CELL OCCLUSIVE

Impede o crescimento epitelial para baixo.

## OPTIMIZED FLEXIBILITY

Rígida o suficiente para a colocação/posicionamento, é facilmente drapeada sobre o cume.



*A construção de fibra reconstituída permite a integração do tecido, evitando a passagem direta de células epiteliais.*